

Clinical Centrifuge

Clinical Centrifuge

A range of centrifuges for general purpose use in Pathology Labs, Small Hospitals, Schools and Colleges. The centrifuge with 5 step speed regulator is provided with Fixed Angle Rotor Head.

Technical Data

Model	Unit	TC 415 A	TC 615 A	TC 815 A	TC 245 A	TC 2415 A
Max. Speed	RPM	4500	4500	4500	4000	4000
Max. RCF	xg	2050	2050	2050	1650	1650
Max. Tube Size	ml	15	15	15	15	1.5/2
Max. Capacity	ml	60	90	120	120	48
Width	mm	210	210	210	210	210
Depth	mm	235	235	235	235	235
Height	mm	230	230	230	255	255
Weight	kg	6	6	6	6	6
Connected Load	kVA	0.11	0.11	0.13	0.13	0.13



TC 415 A