

Magnetic Stirrer (With Hot Plate)

Magnetic Stirrer (With Hot Plate)

This group of Magnetic Stirrers are similar to the above series but have additional hotplate to achieve heating of solutions to be stirred. The heating is controlled by energy regulator which ensures regulation at any desired percentage of the capacity of heater.



MS 203



MS 205

Technical Data

Type	MS 2012	MS 203	MS 205	MS 211	MS 221
Max. speed RPM	1800	1800	1800	1200	1200
Stiring Capacity Ltrs.	1	2	5	10	20
Dimension (WxDxH) mm	140x150x120	160x165x130	190x190x170	250x250x195	340x340x215
Stiring Paddle x L mm	8x25	8x38	9.5x50	9.5x50	13x75
Weight	2.5 Kg.	2.9 Kg.	4.8Kg.	11 Kg.	12 Kg.