

Pathological Centrifuge



TC 725 A

Pathological Centrifuge

Model TC 725 A are centrifuges designed for Pathology Labs, Hospitals, Blood Banks etc. where large number of low volume samples are to be centrifuged. The centrifuge is provided with stepless speed regulator and 1 to 60 minutes preset timer with by pass switch and Fixed Angle Rotor Head to carry 5 ml disposable tubes.

Technical Data		
Model	Unit	TC 725 A
Max. Speed	RPM	4000
Max. RCF	xg	2150
Max. Tube Size	ml	5
Max. Capacity	ml	360
Width	mm	300
Depth	mm	300
Height	mm	250
Weight	kg	10
Connected Load	kVA	0.22