


Vortex Mixer



VM 301

Technical Data
(VM 301)
W: 150 mm x (D): 150 mm x (H): 125 mm
Weight : 3.1 kg, Connected Load : 30 Watts

Vortex Mixer

The unit is meant for quick and thorough mixing of liquids in test tubes, small flasks or bottles. The mixer consists of an enamel finished housing containing the motor with an eccentric drive to a rubber buffer on the top. The tube is held against the vibrating buffer and the intensity of agitation controlled by varying the position in which the tube is held. The thyristorised speed controller enables adjustment of degree of oscillations.

The above unit operates on 220/230 V, 50 cycles, 1 A.C. Supply Dimensions and specifications given above are approximate and not binding. We reserve the right to make improvements modifications without notice.