

Accessories Fibreoptics Led Illuminators




**OP-12R High Intensity
Fluorescence LED illumination
(Front & Back view)**



**OPTOLITE 150
FIBRE OPTIC ILLUMINATOR
WITH FLEXIBLE LIGHT GUIDE
AND POWER SUPPLY**



**EP-3A UNIVERSAL
HALOGEN ILLUMINATOR
with variable light control**



**GS-2 GIPPON-JAPAN
Cross Arm Stand**



**GM-1 GIPPON-JAPAN
Long Arm Stand**

ACCESSORIES FOR

"GIPPON-JAPAN" STEREO AND ZOOM STEREO MICROSCOPE

- GE-100 SPARE EYEPIECE 10X, 15X, 20X (EACH PAIR). CODE: GSE-6060**
GE-101 SPARE EYEPIECE 30X (EACH PAIR). CODE: GSE-6061
GAL-102 SUPPLEMENTARY LENS 0.5X, 0.7X, 1.5X, 2 X (EACH). CODE: SA-6666

FLUORESCENCE RING ILLUMINATORS

- OP-11R FLUORESCENCE RING ILLUMINATOR. CODE: OPR-7061**
GA-101D SPARE CIRCULAR FLUORESCENCE BULB. CODE: SDA-7062
LED RING ILLUMINATOR (CE APPROVED)

OP-12R HIGH INTENSITY LED RING ILLUMINATOR 10W: OP-12R. High Intensity Fluorescent (LED) Ring illuminator, consisting of 64 LED lamps to provide shadow illumination with variable light control. **CODE: OPR-7071**

GA-102L SPARE CIRCULAR RING OF 64 LED BULB 10W (IMPORTED). A Complete Ring with 64 LED Lamps duly mounted for replacement. **CODE: SDA-7064**

"GIPPON" FIBRE OPTIC ILLUMINATORS

(OFFERED WITH LD ILLUMINATORS ON REQUEST)

OP-150S FIBRE OPTIC ILLUMINATOR An elegant halogen cold light illuminator with powerful blower to provide a brilliant spot of cool light, supplied with Single Light Guide with super bright 21V 150W halogen lamp with built in reflector and variable light control. (CE approved). **CODE: OP-9006**

OP-150D FIBRE OPTIC ILLUMINATOR An elegant cold light illuminator, with powerful blower to provide a brilliant spot of cool light, supplied with Double Light Guide with super bright 21V 150W halogen lamp with built in reflector and variable light control. (CE approved). **CODE: OP-9007**

OP-150R FIBRE OPTIC ILLUMINATOR An elegant cold light illuminator, with powerful blower to provide a brilliant spot of cool light, supplied with Ring type Light Guide with super bright 21V 150W halogen lamp with built in reflector and variable light control. (CE approved). **CODE: OP-9008**

SPARE HALOGEN LAMP 21W-150W WITH BUILT IN REFLECTOR.

CODE: OP-9018

EP-3A UNIVERSAL HALOGEN ILLUMINATOR: A sturdy incident halogen illuminator with flexible metal cord for convenient use with all brands and models of Biological and Scientific Microscopes equipped with 6V-20W Halogen Lamp with reflector unit. The illuminator provides bright illumination which is adjustable through variable light control transformer operable on 220V/110V AC. **CODE: EP-9010**

STANDS

GS-2 GIPPON-JAPAN CROSS ARM STAND is an universal cross arm stand with maximum flexibility of movement and can be rotated 360° to examine very large objects directly. Iron soleplate, wearable and stainless steel metal post & rail used with head holding unit and can observe object, enlarge the range of observation. The size of soleplate is 254mm x 254mm. Height of the post 475mm and diameter is 37mm. **CODE: GUS-9009**

GM-1 GIPPON JAPAN LONG ARM STAND is a well designed long arm stand with a platform, specially designed for Electronic industries etc. to observe PCB, Large Circuit Boards and similar objects under the Microscope. The distance between from the center to the main body center is 140mm. Range of adjustment the focus is 100mm (Focus post is 50mm), the diameter of the center working stage is 95mm. The height of the post is 240mm. **CODE: GFM-9016**