

## Penta Opto Multiviewing Microscope



**BIOLUX-NXT-PHC**  
PHASE CONTRAST MICROSCOPE  
offered with Plan Phase objectives  
& Turret Phase Contrast



**TURRET PHASE CONTRAST**



**BIOPLUS-XLS-PHT** Binocular  
Phase Contrast Microscope  
offered with Turret / Slider  
Phase Contrast / Dark Field



**SLIDER PHASE CONTRAST**



**OPTO-DF** Dark Field

### “BIOLUX-NXB/NXT” PHASE CONTRAST MICROSCOPES

**BIOLUX-NXB-PHC UNIVERSAL PHASE CONTRAST CUM BRIGHT FIELD Binocular Microscope with imported sidentopf Binocular Head, having reverse mounted ball bearing quintuple nose piece and co-axial focusing control.** Supplied with DIN PLAN PHASE CONTRAST (TURRET TYPE). Eyepiece: WF-10x (Paired) & Din Plan Phase Objective: PLPH-10x, PLPH20x, PLPH-40x & PLPH-100x & Scanner Objective: PL4x supplied with centering telescope and green filter. **CODE: BIP-2015**

**BIOLUX-NXT-PHC UNIVERSAL PHASE CONTRAST CUM BRIGHT FIELD TRINOCULAR MICROSCOPE.** Same as Biolux-CXB-PHC but supplied with sidentopf Trinocular head (Compensation free). **CODE: BIP-2016**

### “BIOLUX-NXB/NXT” DARK FIELD MICROSCOPES

**BIOLUX-NXB-DF UNIVERSAL DARK FIELD CUM BRIGHT FIELD BINOCULAR MICROSCOPE WITH SIDENTOPF BINOCULAR HEAD AND REVERSE MOUNTED NOSE PIECE.** Supplied with additional Dark Field Condenser. Eyepiece: WF-10x (Paired). Objective: (4 Nos.) SP-4x, SP-10x, SP-40x & SP-100x (With Iris). **CODE: BID-2011**

**BIOLUX-NXT-DF UNIVERSAL DARK FIELD CUM BRIGHT FIELD TRINOCULAR MICROSCOPE.** Same as Biolux-CXB-DF but supplied with sidentopf Trinocular head (Compensation Free). **CODE: BID-2012**

### RESEARCH BINOCULAR MICROSCOPE

#### PHASE CONTRAST MICROSCOPE

**BIOPLUS-XLS-PHS PHASE CONTRAST MICROSCOPE: (SLIDER PHASE 10X AND 40X)** Supplied with Slider Phase Contrast unit, with Semi Plan Phase Objectives SP-PH 10x and SP-PH40x, along with Phase condenser and sliders with annular for 10x and 40x Phase Objective. Standard objective SP4x and SP100x Oil immersion are supplied for normal bright field observation along with paired WF10x eyepiece. Complete with centering telescope and green filter. **CODE: BXL-1775**

**BIOPLUS-XLS-PHT PHASE CONTRAST CUM BRIGHT FIELD BINOCULAR MICROSCOPE (TURRET PHASE WITH FOUR OBJECTIVES)** with high transmission Sidentopf Binocular Head, having reverse mounted quadruple nose piece, imported co-axial focusing control and built in halogen illuminator. Supplied with **DIN STANDARD SEMI PLAN PHASE CONTRAST (TURRET TYPE)** Eyepiece: WF-10x (Paired) & Semi Plan Phase Objective (4 Nos.): SP PH-10x, SP PH20X, SP PH-40x & SP PH-100x. Supplied with centering telescope and green filter. **CODE: BXL-1788**

#### DARK FIELD MICROSCOPE

**BIOPLUS-XL-DF DARK FIELD CUM BRIGHT FIELD BINOCULAR MICROSCOPE WITH HIGH TRANSMISSION BINOCULAR HEAD AND REVERSE MOUNTED NOSE PIECE.** Supplied with interchangeable Dark Field Condenser and Din Standard Semi Plan 100x with IRIS, in place of standard 100x Oil immersion. Eyepiece: WF-10x (Paired). Objective (4 Nos.): SP-4x, SP-10x SP-40x & SP-100x (With Iris). **CODE: BXL1790**

**EXTRA CHARGE FOR SUPPLY WITH TRINOCULAR HEAD IN PLACE OF BINOCULAR HEAD IN ALL MODELS OF BIOLUX SERIES MICROSCOPES**