

INCUBATOR



Construction : Double walled

Inside SS / GI

Outside SS / GI

Controls : Thermostatic / Digital

Accuracy : $\pm 1^{\circ}\text{C}$ / $\pm 0.5^{\circ}\text{C}$

Temp. Range : 5°C > ambient to Max 60°C

Size :	Vol	Weight
12 X 12 X 12"	27 ltrs	28 Kgs
14 X 14 X 14"	42 ltrs	33 Kgs
18 X 18 X 18"	91 ltrs	51 Kgs
18 X 18 X 24"	121 ltrs	66 Kgs
24 X 24 X 24"	216 ltrs	77 Kgs
24 X 24 X 36"	324 ltrs	85 Kgs