


## Magnetic Stirrer (without Heating)



**MS 202**

### Magnetic Stirrer (Without Heating)

These equipments are intended for stirring solutions of low viscosity in open or closed vessels. Stirring is accomplished by a magnetized teflon coated paddle to be placed in the container. Each apparatus of this group can be operated on a bench, a stand or a laboratory scaffold as well. The intensity of stirring, depending on requirements, can be adjusted with the thyristorised voltage controller incorporated which also ensures step less adjustments from very low to the maximum agitation.

**Technical Data**

Type	MS 202	MS 204	MS 210	MS 220	MS 250
Max. speed RPM	1800	1800	1200	1200	1200
Stirring Capacity Ltrs.	2	5	10	20	50/75
Dimension (WxDxH)	160x165x120	190x190x170	250x250x170	340x340x195	450x450x250
Stirring Paddle x L mm	8x38	9.5x50	9.5x50	13x75	13x100
Weight	1.5 Kg.	3.8 Kg.	10 Kg.	11 Kg.	25 Kg.

MAGNETIC STIRRERS up to 200 litres capacity can be manufactured and supplied to suit users specific requirements.