

KF KFS Medical Microscope

MONOCULAR MICROSCOPE

KF-33 GIPPON-GETNER ADVANCE STUDENT MICROSCOPE A study Monocular Microscopes manufactured in collaboration with Gippon-Japan, & equipped with Abbe condenser and iris diaphragm moveable on Rack & Pinion. Supplied with all imported color coded antifungus "Gippon-Getner" optics including Eyepiece: (2 No) WF10x and H5x/ H15x ACHROMAT OBJECTIVE:10x & 45x (Spring Mount) (Thermocole Packing).

CODE: KF 1033

KF-4 (i) GIPPON-GETNER MEDICAL MICROSCOPE WITH BUILT IN CO-AXIAL MECHANICAL STAGE with abbe condenser 1.25 N.A. moveable on Rack & Pinion with iris diaphragm, filter holder and built in co-axial ball bearing graduated mechanical stage with low position controls. Imported, Antifungus color coded Gippon-Getner, opics including Eyepiece: (1 No) WF 10x. Objectives : (3 Nos.) 10x, 40 / 45x and 100x oil immersion (spring mount), (Thermocole Packing).

CODE: KF-1004

ACCESSORIES

GA-405 SPARE "GIPPON-GETNER" H5x or H15x Eyepiece for KF Series Microscopes. (Any 1 No.)

CODE: KA-1030

GA-501 SPARE "GIPPON-GETNER" 4x objective for KF & BIOPLUS Microscopes

CODE: KA-1031

GL-1 SUB STAGE MICROSCOPE LAMP Now upgraded and supplied with imported compact PIGMY LAMP for centered illumination.

CODE: KA-1050

WF 105 SPARE LAMP 15W PIGMY LAMP (IMPORTED) for GL-1 Sub Stage lamp

CODE: KA-1051

WF-101 PLYWOOD STORING CABINET.

CODE: KA-1040

INCLINED MONOCULAR & PHOTO MICROSCOPES

BIOPLUS-M INCLINED MONOCULAR RESEARCH MICROSCOPE, WITH HIGH TRANSMISSION MONOCULAR TUBE, REVERSE MOUNTED QUADRUPLE NOSE PIECE and "GIPPON-GETNER" ACHROMAT optics (Anti fungus). Eyepiece WF-10x Objectives: (4 Nos.) 4x, 10x, 40x & 100x (Oil Immersion).

CODE: BXL-1665

BIOLINE-MP INCLINED MONOCULAR PHOTO MICROSCOPE - A unique model of Monocular Microscope equipped with built in phototube for photomicrography and use of video equipments. Highly economical model, most competitively priced and yet providing all facilities for various Microscopy applications and Cameras and CCD Cameras etc. Supplied with Anti-Fungus hard coated "GIPPON-GETNER" Din Achromat Optics (Anti fungus). Eyepiece Wf-10x (Paired) Objectives (4 Nos.) 4x, 10x, 40x & 100x (Oil immersion).

CODE: BXL-1777

BINOCULAR/TRINOCULAR MICROSCOPES

BIOPLUS-XL BINOCULAR RESEARCH MICROSCOPE, WITH HIGH TRANSMISSION BINOCULAR HEAD, REVERSE MOUNTED QUADRUPLE NOSE PIECE and "GIPPON-GETNER" ACHROMAT OPTICS (Anti fungus). Eyepiece WF-10x (Paired), Objectives (4 Nos.): 4x, 10x, 40x & 100x (Oil Immersion).

CODE: BXL-1775

BIOPLUS-XLT TRINOCULAR RESEARCH MICROSCOPE, WITH HIGH TRANSMISSION TRINOCULAR HEAD, REVERSE MOUNTED QUADRUPLE NOSE PIECE and "GIPPON-GETNER" ACHROMAT OPTICS (Anti fungus). Eyepiece WF-10x (Paired), Objectives (4 Nos.): 4x, 10x, 40x & 100x (Oil Immersion).

CODE: BXL-1776



KF-4 MEDICAL MICROSCOPE
(Ball Bearing Co-axial Stage)



GL-1 SUBSTAGE LAMP



BIOLINE-M
MONOCULAR MICROSCOPE



BIOLINE-MP
MONOCULAR MICROSCOPE



BIOPLUS-XLT
Trinocular Microscope
With Sidentopf Head Semi Plan Optics.