



OVEN



Construction : Double walled

Inside SS / GI

Outside SS / GI

Controls : Thermostatic / Digital

Accuracy : $\pm 1^\circ\text{C}$ / $\pm 0.5^\circ\text{C}$

Temp. Range : 5°C > ambient to Max 250°C

Size :	Vol	Weight
12 X 12 X 12"	27 ltrs	30 Kgs
14 X 14 X 14"	42 ltrs	35 Kgs
18 X 18 X 18"	91 ltrs	55 Kgs
18 X 18 X 24"	121 ltrs	70 Kgs
24 X 24 X 24"	216 ltrs	80 Kgs
24 X 24 X 36"	324 ltrs	90 Kgs



LAB ENTERPRISE

30/237, dadabhai road, vile parle (w), mumbai - 400056.
contact no - 9769803350, telefax - 2671 5001
e-mail - info@labenterprise.com, labentp@gmail.com.

