

VACUUM OVEN ROUND

MODEL EHT - 15 (A)



Round Model

Round Models

22.5 cms Dia X 30.0 cms Depth
30.0 cms Dia X 37.5 cms Depth
37.5 cms Dia X 45.0 cms Depth
45.0 cms Dia X 60.0 cms Depth

Double Walled, outer made of M.S. powder coated & inner wall made of stainless steel in cylindrical / Rectangular shape. Stainless steel suspension type door closing mechanism with glass windows is provided. Heaters are clamped on to the outside of the body controlled by Thermostat / Digital Display range 50 °C to 200 °C

A vacuum type stop cock and a vacuum gauge are fitted on the top. Top shelves are provided. (L shaped thermometer fix on shelves in case of Thermostatic Model only)